

آسان تجهیز ایرانیان
Asan Tajhiz Iranian



ست ابزار کایفوپلاستی

Kyphoplasty Smaller Instrument System



www.atico-trade.com

شرکت آسان تجهیز ایرانیان
نماینده رسمی و انحصاری محصولات
کایفوپلاستی کمپانی BMK در ایران



- GUARDIAN-SG™-Kyphoplasty Smaller Instrument System
- BARCLAY®-Radiopaque Bone cement



GUARDIAN-SG™ KYPHOPLASTY SMALLER INSTRUMENT SYSTEM



ست ابزار کایفوپلاستی Guardian-SG، شامل ابزار ویژه ستون مهره‌ها می‌باشد که برای دسترسی آسانتر و جلوگیری از سرعت پیشرفت شکستگی ستون مهره‌ها و ایجاد فضای خالی (جهت ترمیم و پر نمودن شکستگی با تزریق سیمان) در بافت استخوانی اسفنجی طراحی شده است. این سیستم، ویژگی‌های مذکور را از طریق ابزار با قطر کوچک (Small Gauge) ارائه می‌دهد.

ویژگی خاص ابزار کایفوپلاستی Guardian-SG

ابزار کانولا (Working Cannula) جهت دسترسی سریع و کارآمد از طریق پوست به استخوان می‌باشد. این شرکت، کانولا را با قطر کوچک (10G) برای کاهش تروما و درد طراحی نموده است که موجب کاهش صدمات به بیمار حین عمل می‌گردد.



Beveled Tip



Diamond tip



※ Bone biopsy device



※ Working Cannula of 10gauge

Component of GUARDIAN-SG TM



• Accessory

- Balloon Catheter(T)	- Drill(T)
- Compressor	- Cement Infusion(T)
- Cannula(T)	- Cement Bar(T)
- Trocar-S(T)	

GUARDIAN-SG TM Ordering Information

► Direct access Version

Balloon No.	Article No.	Balloon No.	Article No.
One Balloon	BMK(T)10αS	Two Balloon	2BMK(T)10αS
	BMK(T)15 αS		2BMK(T)15 αS
	BMK(T)20 αS		2BMK(T)20 αS

Balloon Catheter Description

Article No.	Length	PSI Rating	Volume(cc)
BMK(T)10α	10	400	4
BMK(T)15 α	15	400	4
BMK(T)20 α	20	400	6

*This product is sterilized by EO gas and followed CE standard regulations.

Barclay[®]

Radiopaque Bone cement

Barclay[®] is a self-hardening, ready to use, high-viscosity bone cement. Barclay[®] contains a high percentage of radiopaque agent(%30) for maximum visibility and radiological control during a percutaneous kyphoplasty procedure.



Features

- High radiopacity for maximum visibility and radiological control
- High viscosity enables controlled cement delivery
- Longer working time, allows for safe injection during kyphoplasty
- Lower degree of exothermic reaction

Ordering Information

Article No.	Description
KB20	20g powder + 9.2g Liquid
KB25	25g powder + 9.7g Liquid
KB30	30g powder + 13.8g Liquid

*KB 21: "MUDLOCK" brand name

► Handling Times for Barclay Bone cement at °23c

